

2-Amino-3-pyridinecarbonitrile

Art. ID	TRC-A629155-5G
Unit	5 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2811 III E1

Description

Appearance: Light Beige to Dark Beige Solid /// Application Notes: 2-Amino-3-pyridinecarbonitrile is used as a reagent in the synthesis of novel imidazo[1,2-a]pyridine derivatives as potent c-Met inhibitors and antimycobacterials. /// References: Li, C., et al.: ACS Med. Chem. Lett., 6, 507 (2015)| Ramachandran, S., et al.: Bioorg. Med. Chem. Lett., 23, 4996 (2013)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-3-pyridinecarbo nitrile	[24517-64-4]				