

(3-Aminopropyl)trimethoxysilane

Art. ID	TRC-A629320-25ML
Unit	25 ml
Deliverydetails	No Dangerous Good /not restricted

Description

Appearance: Clear Colourless Oil /// Application Notes: (3-Aminopropyl)trimethoxysilane is used in the synthesis of aspirin imprinted polymer on silica gel surface and appreciation on drug sustained release. Fluorine-doped, lectin functionalized ZnO nanoarrays as 3D nano-biointerface for bacterial detection. /// References: Lu, P., et al.: Gaofenzi Cailiao Kexue Yu Gongcheng, 33, 130 (2017)| Zheng, L., Talanta, 167, 600 (2017)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(3-Aminopropyl)trimethoxysilane	[13822-56-5]				