

3-Aminopiperidine-2,6-dione Hydrochloride

Art. ID TRC-A632980-5G
Unit 5 g
Deliverydetails No Dangerous Good

Description

Appearance: Off-White Solid /// Application Notes: 3-Aminopiperidine-2,6-dione hydrochloride is a reagent for preparing lenalidomide that can induce ubiquitination and degradation of CK1alpha in del(5q) MDS. It can also be used to prepare phthalimide conjugates that can promote ligand-dependent target protein degradation. /// References: Kronke, J., et al.: Nature 523, 183 (2015)| Liu C., et al.: Faming Zhuanli Shenqing CN 105348257 A 20160224 (2016)| Rao, R., et al. : PCT Int. Appl. WO 2016024286 A2 20160218 (2016)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Aminopiperidine-2,6-dione Hydrochloride	[24666-56-6]				