

Ammonium Dichromate

Art. ID TRC-A633980-25G
Unit 25 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1439 II E2

Description

Application Notes: Ammonium Dichromate is used in pyrotechnics as well as in the manufacture of photosensitive films used in electronical devices. Ammonium dichromate acting as the photoactive site. /// References: Harvard, J. et al.: Chem Mater. 11, 719 (1999)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium dichromate	[7789-09-5]				