

Ampicillin Trimer Trisodium Salt

Art. ID	TRC-A634395-2.5MG
Unit	2,5 mg
Deliverydetails	Dry ice shipment

Description

Appearance: White Solid /// Application Notes: Ampicillin Trimer Trisodium Salt can be used to develop HPLC methods that detect compounds or impurities related to ampicillin in bulk drug. Ampicillin Trimer Trisodium Salt is an impurity of Ampicillin (A634300), a beta-lactam antibiotic used for treating bacterial infection. Ampicillin interferes with cell wall synthesis in bacteria. /// References: Raju, Ch. B. V. N., et al.: Anal Chem : Ind J, 7, 354 (2008)| Fisher, J. F., et al.: Chem. Rev. 105, 395 (2005)| Sill, M. L., et al.: Antimicrob. Agents Chemother. 52, 1551 (2008)| Ivashkiv, E., et al.: Anal. Profiles Drug Subs. 2, 1 (1973)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ampicillin Trimer Triso dium Salt					