

## **p-Anisidine**

Art. ID TRC-A673505-50G  
Unit 50 g  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2811 III E1

### **Description**

Category: Amines, Aromatics, Miscellaneous Reagents /// Appearance: Dark Grey to Black Solid /// Application Notes: p-Anisidine, a grey-brown solid, is the most toxic of the three isomers of anisidine. It can be used for the detection of fats and oils, since p-Anisidine reacts with aldehydes and ketones in fats and oils to form products that absorb at 350 nm wavelength of light. /// References: Weast, R.C., CRC. Handbook Chem. Phys. (1981)| Steele, R., Understanding and Measuring the Shelf-Life of Food (2004)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Anisidine	[104-94-9]				