

Aprotinin Hydrochloride

Art. ID TRC-A735430-250MG
Unit 250 mg
Deliverydetails No Dangerous Good

Description

Category: Amino Acids and Derivatives, Peptides /// Appearance: Off-White Solid /// Application Notes: Aprotinin is used in cardiac surgery to minimize bleeding and reduce exposure to blood products. /// References: Mangano, D. et al.: New Eng. J. Med., 354, 353 (2006)| Fergusson, D. et al.: New. Eng. J. Med., 358, 2319 (2008)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aprotinin Hydrochloride					