

N-Arachidonoyltaurine

Art. ID TRC-A765080-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Category: Amines, Sulfur and Selenium Compounds /// Application Notes: N-Arachidonoyltaurine is a polyunsaturated lipid that inhibits the Cav3.3 current in T-type calcium channels. /// References: Cazade, M., et. al.: Mol. Pharmacol., 85, 218 (2014)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Arachidonoyltaurine	[119959-65-8]				