

Arsenic Oxide (As₂O₃)

Art. ID TRC-A774130-25G

Unit 25 g

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1561 II E4

Description

Category: Miscellaneous Reagents /// Appearance: White to Pale Yellow Solid /// Application Notes: Arsenic Oxide ((As₂O₃)) is used as a precursor to forestry products, in colorless glass production, and in electronics. In combination with copper(II), acetate arsenic trioxide gives Copper Acetoarsenite (C685255) as a rodenticide. /// References: Grund, S., et al.: Ullmann's Encyclopedia of Industrial Chemistry, Weinheim: Wiley-VCH (2005)| Gaines, T.: Toxicol. Appl. Pharm., 2, 88 (1960)| Jameson, F.: Paint Manu., 13, 6 (1943)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (III) oxide	[1327-53-3]				