

L-Aspartic Acid-d3

Art. ID TRC-A790023-50MG
Unit 50 mg
Deliverydetails No Dangerous Good

Description

Category: Amino Acids and Derivatives, Isotope Labelled Compounds, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White to Pale Beige Solid /// Application Notes: Labelled analogue of L-Aspartic Acid, a non-essential amino acid found in food sources and dietary supplements. L-Aspartic Acid is one of the 20 proteinogenic amino acids, the building blocks of proteins. Its conjugate base L-aspartate is an excitatory neurotransmitter in the central nervous system. /// References: Ager, D.J. et al.: J. Mol. Catal. Enzym., 11, 199 (2001)| Velisek, J. et al.: Cze. J. Food. Sci., 24, 1 (2006)| Nagahama, M. et al.: Cell Tiss. Res., 318, 483 (2004)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Aspartic Acid-d3	[3842-25-9]				