

L-Aspartic Acid Dibenzyl Ester p-Toluenesulfonate Salt

Art. ID	TRC-A790070-10G
Unit	10 g
Deliverydetails	No Dangerous Good

Description

Category: Amines, Aromatics, Chiral Reagents, Miscellaneous Reagents /// Appearance: White Solid /// Application Notes: L-Aspartic Acid Dibenzyl Ester p-Toluenesulfonate Salt is used in the preparation of Betulinic Acid derivatives that show inhibition towards HIV. It is also used in the preparation of piperidine functional amides which act as potent and highly selective human beta3 agonists. /// References: Soler, F. J. Med. Chem., 39, 1069 (1999)]

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Aspartic Acid Dibenzyl Ester p-Toluenesulfonate Salt	[2886-33-1]				