

Auramine O Hydrochloride (~75%)

Art. ID TRC-A794545-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2811 III E1

Description

Appearance: Dark Yellow to Orange Solid /// Application Notes: Auramine O is a substituted diphenyl-methane cationic dye. Auramine O has been used as a probe for studies of microheterogeneous systems, from micelles, proteins, polymers to silica surfaces. /// References: Ferreira, A.U.C., et al.: J. Lumin., 136, 63 (2012)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Auramine O	[2465-27-2]				