

## **Acetonitrile**

Art. ID TRC-A163890-100ML  
Unit 100 ml  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1648 II E2

### Description

Category: Intermediates, Miscellaneous Reagents /// Appearance: Clear Colourless Liquid /// Application Notes: Acetonitrile is used as a polar aprotic solvent in organic synthesis. /// References: Hu, H.H., et al.: Phy. Chem. Phys., 15, 7666 (2013)| Kim, J.Y., et al.: J. Mater. Chem. A., et al.: 1, 5982 (2013)|

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | Acetonitrile      | [75-05-8]  |                     |      |        |        |