

Acetophenone-d3

Art. ID	TRC-A164023-500MG
Unit	500 mg
Deliverydetails	No Dangerous Good /not restricted

Description

Application Notes: Acetophenone-d3 is an isotopic analog of acetophenone (A164015), a reagent used in the production of fragrances and resin polymers. This compound is a contaminant of emerging concern (CECs) //
References: Perez, D., et al. | Bioorg. Med. Chem., 17, 6914 (2009), Giet, J., et al.: J. Agric. Food Chem., 57, 6787 (2009), Dimoula, K., et al.: Biotechnol. J., 4, 712 (2009)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetophenone-beta,beta, beta-d3	[17537-31-4]				