

Acetic Acid-d4

Art. ID TRC-A167657-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2789 3 II E2

Description

Category: Isotope Labelled Compounds, Labeling and Diagnostics Reagents /// Appearance: Clear Colourless Liquid /// Application Notes: Labelled Acetic Acid is used for analysis and profiling of biological samples resolved through NMR spectroscopy. It helps with the labelling of RNA for high resolution NMR. /// References: Cahoreau, E. et al.: Anal. Bio., 427, 158 (2012)|

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Acetic acid-d4 | [1186-52-3] | | | | |