

N-Acetyl Lorcaserin-d3

Art. ID TRC-A185202-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Application Notes: N-Acetyl Lorcaserin-d3 is the isotope labelled analog of N-Acetyl Lorcaserin (A185200), an impurity of Lorcaserin (L469890, HCl) which is a novel selective 5-HT_{2C}-receptor agonist for the treatment of obesity. N-Acetyl Lorcaserin is also used as a reagent in the synthesis of isoquinoline, isoindoline, and benzoazepine amides as therapeutic modulators of the histamine H₃ receptor. /// References: Pauli, M., et al.: Formulary, 45, 180 (2010)| Ioannides-Demos, L.L., et al.: J. Obesity, 18, 179674 (2011)| Santora, V.J., et al.: PCT Int. Appl. WO 2009105206 A1 20090827. Aug 27, 2009

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl Lorcaserin-d3					