

## **CONOSTAN - Multi-element standard in 75 blank oil: S21 + K+Li: 815 ppm**

Art. ID                      CON-150-021-065  
Unit                         200 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0815	%		
	Aluminium (Al)	[7429-90-5]	815	µg/g		
	Barium (Ba)	[7440-39-3]	0,0815	%		
	Barium (Ba)	[7440-39-3]	815	µg/g		
	Boron (B)	[7440-42-8]	0,0815	%		
	Boron (B)	[7440-42-8]	815	µg/g		
	Cadmium (Cd)	[7440-43-9]	0,0815	%		
	Cadmium (Cd)	[7440-43-9]	815	µg/g		
	Calcium (Ca)	[7440-70-2]	0,0815	%		
	Calcium (Ca)	[7440-70-2]	815	µg/g		
	Chromium (Cr)	[7440-47-3]	0,0815	%		
	Chromium (Cr)	[7440-47-3]	815	µg/g		
	Copper (Cu)	[7440-50-8]	0,0815	%		
	Copper (Cu)	[7440-50-8]	815	µg/g		
	Iron (Fe)	[7439-89-6]	0,0815	%		
	Iron (Fe)	[7439-89-6]	815	µg/g		
	Lead (Pb)	[7439-92-1]	0,0815	%		
	Lead (Pb)	[7439-92-1]	815	µg/g		
	Lithium (Li)	[7439-93-2]	0,0815	%		
	Lithium (Li)	[7439-93-2]	815	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,0815	%		
	Magnesium (Mg)	[7439-95-4]	815	µg/g		
	Manganese (Mn)	[7439-96-5]	0,0815	%		
	Manganese (Mn)	[7439-96-5]	815	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,0815	%		
	Molybdenum (Mo)	[7439-98-7]	815	µg/g		
	Nickel (Ni)	[7440-02-0]	0,0815	%		
	Nickel (Ni)	[7440-02-0]	815	µg/g		
	Phosphorus (P)	[7723-14-0]	0,0815	%		
	Phosphorus (P)	[7723-14-0]	815	µg/g		
	Potassium (K)	[7440-09-7]	0,0815	%		

Potassium (K)	[7440-09-7]	815	µg/g
Silicon (Si)	[7440-21-3]	0,0815	%
Silicon (Si)	[7440-21-3]	815	µg/g
Silver (Ag)	[7440-22-4]	0,0815	%
Silver (Ag)	[7440-22-4]	815	µg/g
Sodium (Na)	[7440-23-5]	0,0815	%
Sodium (Na)	[7440-23-5]	815	µg/g
Tin (Sn)	[7440-31-5]	0,0815	%
Tin (Sn)	[7440-31-5]	815	µg/g
Titanium (Ti)	[7440-32-6]	0,0815	%
Titanium (Ti)	[7440-32-6]	815	µg/g
Vanadium (V)	[7440-62-2]	0,0815	%
Vanadium (V)	[7440-62-2]	815	µg/g
Zinc (Zn)	[7440-66-6]	0,0815	%
Zinc (Zn)	[7440-66-6]	815	µg/g