

CONOSTAN - Multi-element standard in 75 blank oil: S21 + K+Sb: 30 ppm

Art. ID CON-150-021-076
Unit 200 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,003 | % | | |
| | Aluminium (Al) | [7429-90-5] | 30 | µg/g | | |
| | Antimony (Sb) | [7440-36-0] | 0,003 | % | | |
| | Antimony (Sb) | [7440-36-0] | 30 | µg/g | | |
| | Barium (Ba) | [7440-39-3] | 0,003 | % | | |
| | Barium (Ba) | [7440-39-3] | 30 | µg/g | | |
| | Boron (B) | [7440-42-8] | 0,003 | % | | |
| | Boron (B) | [7440-42-8] | 30 | µg/g | | |
| | Cadmium (Cd) | [7440-43-9] | 0,003 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 30 | µg/g | | |
| | Calcium (Ca) | [7440-70-2] | 0,003 | % | | |
| | Calcium (Ca) | [7440-70-2] | 30 | µg/g | | |
| | Chromium (Cr) | [7440-47-3] | 0,003 | % | | |
| | Chromium (Cr) | [7440-47-3] | 30 | µg/g | | |
| | Copper (Cu) | [7440-50-8] | 0,003 | % | | |
| | Copper (Cu) | [7440-50-8] | 30 | µg/g | | |
| | Iron (Fe) | [7439-89-6] | 0,003 | % | | |
| | Iron (Fe) | [7439-89-6] | 30 | µg/g | | |
| | Lead (Pb) | [7439-92-1] | 0,003 | % | | |
| | Lead (Pb) | [7439-92-1] | 30 | µg/g | | |
| | Magnesium (Mg) | [7439-95-4] | 0,003 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 30 | µg/g | | |
| | Manganese (Mn) | [7439-96-5] | 0,003 | % | | |
| | Manganese (Mn) | [7439-96-5] | 30 | µg/g | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,003 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 30 | µg/g | | |
| | Nickel (Ni) | [7440-02-0] | 0,003 | % | | |
| | Nickel (Ni) | [7440-02-0] | 30 | µg/g | | |
| | Phosphorus (P) | [7723-14-0] | 0,003 | % | | |
| | Phosphorus (P) | [7723-14-0] | 30 | µg/g | | |
| | Potassium (K) | [7440-09-7] | 0,003 | % | | |

| | | | |
|---------------|-------------|-------|------|
| Potassium (K) | [7440-09-7] | 30 | µg/g |
| Silicon (Si) | [7440-21-3] | 0,003 | % |
| Silicon (Si) | [7440-21-3] | 30 | µg/g |
| Silver (Ag) | [7440-22-4] | 0,003 | % |
| Silver (Ag) | [7440-22-4] | 30 | µg/g |
| Sodium (Na) | [7440-23-5] | 0,003 | % |
| Sodium (Na) | [7440-23-5] | 30 | µg/g |
| Tin (Sn) | [7440-31-5] | 0,003 | % |
| Tin (Sn) | [7440-31-5] | 30 | µg/g |
| Titanium (Ti) | [7440-32-6] | 0,003 | % |
| Titanium (Ti) | [7440-32-6] | 30 | µg/g |
| Vanadium (V) | [7440-62-2] | 0,003 | % |
| Vanadium (V) | [7440-62-2] | 30 | µg/g |
| Zinc (Zn) | [7440-66-6] | 0,003 | % |
| Zinc (Zn) | [7440-66-6] | 30 | µg/g |