

## **Uric acid**

Art. ID	NIST-913b
Unit	10 g
Deliverydetails	No Dangerous Good /not restricted

### **Description**

This Standard Reference Material (SRM®) is certified as a chemical of known purity. It is intended primarily for use in the calibration and standardization of procedures for uric acid determinations employed in clinical analysis and for routine critical evaluation of the daily working standards used in these procedures. A unit of NIST-913b consists of one bottle containing 10 g of crystalline uric acid. Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified purity (mass fraction)	Uric Acid	[69-93-2]	99,8 ± 0,2%			