

Tin-base Lead-free Solder, Cast, Certified Reference Material (CRM), disc 40 mm x 15 mm

Art. ID MBH-74XAMF
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,001	%		
	Silver (Ag)	[7440-22-4]	0,496	%		
	Arsenic (As)	[7440-38-2]	0,0038	%		
	Bismuth (Bi)	[7440-69-9]	0,190	%		
	Cadmium (Cd)	[7440-43-9]	0,0061	%		
	Copper (Cu)	[7440-50-8]	3,07	%		
	Indium (In)	[7440-74-6]	0,0082	%		
	Iron (Fe)	[7439-89-6]	~0,001	%		
	Nickel (Ni)	[7440-02-0]	0,0260	%		
	Phosphorus (P)	[7723-14-0]	~0,002	%		
	Lead (Pb)	[7439-92-1]	0,126	%		
	Antimony (Sb)	[7440-36-0]	1,064	%		
	Selenium (Se)	[7782-49-2]	~0,001	%		