

PAK - Mix 21, each 100 µg/ml in Toluene

Art. ID LM24-SL-10013
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene	[83-32-9]	100	µg/mL		
	Acenaphthylene	[208-96-8]	100	µg/mL		
	Anthracene	[120-12-7]	100	µg/mL		
	Benzo(a)anthracene	[56-55-3]	100	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	100	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	100	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	100	µg/mL		
	Benzo(a)pyrene	[50-32-8]	100	µg/mL		
	Chrysene	[218-01-9]	100	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	100	µg/mL		
	Fluoranthene	[206-44-0]	100	µg/mL		
	Fluorene	[86-73-7]	100	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	100	µg/mL		
	Naphthalene	[91-20-3]	100	µg/mL		
	Phenanthrene	[85-01-8]	100	µg/mL		
	Pyrene	[129-00-0]	100	µg/mL		
	Benzo(e)pyrene	[192-97-2]	100	µg/mL		
	Benzo(j)fluoranthene	[205-82-3]	100	µg/mL		
	Naphthalene D8	[1146-65-2]	100	µg/mL		
	Pyrene D10	[1718-52-1]	100	µg/mL		
	Perylene D12	[1520-96-3]	100	µg/mL		