

PAK - Mix 3, each 100 µg/mL in Toluene

Art. ID LM24-CUS-10015
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Naphthalene D8	[1146-65-2]	100	µg/mL		
	Pyrene D10	[1718-52-1]	100	µg/mL		
	Perylene D12	[1520-96-3]	100	µg/mL		