

**PAK - Mix 3, each 100 µg/mL in Toluene**

Art. ID LM24-CUS-10015  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Naphthalene D8    | [1146-65-2] | 100                 | µg/mL |        |        |
|                  | Pyrene D10        | [1718-52-1] | 100                 | µg/mL |        |        |
|                  | Perylene D12      | [1520-96-3] | 100                 | µg/mL |        |        |