

Tetranic acid, C4:0, Purity: >99%

Art. ID LA-10-0400-17
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Synonyms: Butyric acid / Butanoic acid

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Tetranic acid | [107-92-6] | | | | |