

Cotton Bayer CropScience Non-Modified Leaf DNA

Art. ID AOCS-0306-A3
Unit approx. 10 µg
Deliverydetails No Dangerous Good /not restricted

Description

AOCS-0306-A3 has been prepared by AOCS from cotton genomic DNA extracted from clean leaves. AOCS-0306-A3 is available in sterile, 0.5ml skirted screw-cap self-sealing tubes containing approximately 10 µg of cotton genomic DNA extracted from clean leaves. This reference material has been produced from genomic DNA extracted from clean leaves of non-modified cotton delivered by Bayer BioScience NV. The non-modified cotton used in the preparation of AOCS ISO Guide 34 A2LA Certificate 3438.01 Certificate of Analysis for 0306-A3 Page 2 of 4, © AOCS, 2013 0306-A3 resulted from several cycles of self-pollination of the conventional line.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Presence: No	LLCotton25, GHB614					
Pure, non-modified cotton						