

Copper base, Type: Bronze - Aluminum, XRF Setting-up-sample, disc 40 x 15 mm - please ask for availability

Art. ID SUS-RC33XRF
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	<0,001	%		
	Aluminium (Al)	[7429-90-5]	9,39	%		
	Arsenic (As)	[7440-38-2]	0,02	%		
	Bismuth (Bi)	[7440-69-9]	0,005	%		
	Cadmium (Cd)	[7440-43-9]	0,003	%		
	Cobalt (Co)	[7440-48-4]	0,01	%		
	Chromium (Cr)	[7440-47-3]	0,04	%		
	Copper (Cu)	[7440-50-8]	81,43	%		
	Iron (Fe)	[7439-89-6]	4,2	%		
	Magnesium (Mg)	[7439-95-4]	<0,001	%		
	Manganese (Mn)	[7439-96-5]	0,2	%		
	Nickel (Ni)	[7440-02-0]	4,5	%		
	Phosphorus (P)	[7723-14-0]	<0,01	%		
	Lead (Pb)	[7439-92-1]	<0,01	%		
	Sulfur (S)	[7704-34-9]	0,001	%		
	Antimony (Sb)	[7440-36-0]	<0,001	%		
	Silicon (Si)	[7440-21-3]	0,04	%		
	Tin (Sn)	[7440-31-5]	<0,01	%		
	Titanium (Ti)	[7440-32-6]	<0,001	%		
	Zinc (Zn)	[7440-66-6]	0,09	%		