

10(Z),13(Z),16(Z)-Docosatrienoic acid, C22:3, Purity: >98%, solution in Ethanol, concentration: 10 mg/mL

Art. ID LA-10-2213-4
Unit 5 mg
Deliverydetails De Minimis UN1170 III E1

Description

Synonyms: cis-10,13,16-Docosatrienoic acid / 10-cis,13-cis,16-cis-Docosatrienoic acid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	10(Z),13(Z),16(Z)-Docosatrienoic acid	[18104-42-2]	10	mg/mL		