

4(Z),10(Z),13(Z),16(Z)-Docosatetraenoic acid, C22:4, Purity: >99%, solution in Ethanol, concentration: 50 mg/mL

Art. ID LA-10-2214-5
Unit 10 mg
Deliverydetails De Minimis UN1170 III E1

Description

Synonyms: cis-4,10,13,16-Docosatetraenoic acid / 4-cis,10-cis,13-cis,16-cis-Docosatetraenoic acid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4(Z),10(Z),13(Z),16(Z)- Docosatetraenoic acid	[122068-08-0]	50	mg/mL		