

Custom Solution - 6 components in Methanol

Art. ID LM24-CUS-70179
Unit 5 mL (ampoule)
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Chlorophenol	[95-57-8]	300	µg/mL		
	3-Chlorophenol	[108-43-0]	300	µg/mL		
	4-Chlorophenol	[106-48-9]	300	µg/mL		
	2,4,5-Trichlorophenol	[95-95-4]	30	µg/mL		
	3,4,5-Trichlorophenol	[609-19-8]	30	µg/mL		
	Pentachlorophenol	[87-86-5]	10	µg/mL		