

### **N-Acetyl-4-aminosalicylic Acid-d3**

Art. ID TRC-A168327-50MG  
Unit 50 mg  
Deliverydetails No Dangerous Good

#### Description

Application Notes: N-Acetyl-4-aminosalicylic Acid-d3 is the deuterated form of N-Acetyl-4-aminosalicylic Acid (A168325), which is a versatile reactant used in the preparation of salicylanilide derivatives as inhibitors of the protein tyrosine kinase epidermal growth factor receptor. /// References: Liechti, C., et al.: Eur. J. Med. Chem., 39, 11 (2004)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-4-aminosalicylic Acid-d3					