

## **Acetylenedicarboxylic Acid**

Art. ID	TRC-A173870-50G
Unit	50 g
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2811 III E1

### Description

Application Notes: Acetylenedicarboxylic Acid is used as an electrolyte in the fabrication of anodic nanoporous alumina. It is also used as a reagent in the synthesis of several organic compounds including that of derivatives of metronodazole (antibiotic and antiprotozoal medication) which display increased antimicrobial and antiparasitic properties. /// References: Kikuchi, T., et al.: Electrochem. Lett., 3, C25 (2014)| Al-Masoudi, N., et al.: Monatsh. Chem., 147, 383 (2016)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetylenedicarboxylic a cid	[142-45-0]				