

N-Acetyl-L(-)-leucine

Art. ID TRC-A185035-25G
Unit 25 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Miscellaneous Reagents, /// Appearance: White Solid /// Application Notes: N-Acetyl-L(-)-leucine is used in the preparation of small molecule inhibitors of anti-apoptotic Bcl-2 family proteins. Also used in the preparation of amphiphilic copolymers involving hydrophobic amino acid and oligopeptide side chains for optical tumor imaging in vivo. /// References: Zheng, C. et al.: Bioorg. Med. Chem. Lett., 22, 39 (2012)| Miki, K. et al.: Bull. Chem. Soc. Jap., 85, 1277 (2012)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-L-Leucine - US P	[1188-21-2]				