

o-Acetyltoluene

Art. ID TRC-A189000-100G
Unit 100 g

Description

Category: Aromatics, Miscellaneous Reagents /// Appearance: Clear Colourless Oil /// Application Notes:
o-Acetyltoluene can be used as a reagent to prepare diaryl dihydropyrazoles as inhibitors of mitotic kinesin spindle protein (1). It is also used to synthesize aryl- and indolyl-cyanoaminopyridines with anti-inflammatory, analgesic, and antipyretic activities (2).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | o-Acetyltoluene | [577-16-2] | | | | |