

Adenosine 5'-Triphosphate Disodium Salt Hydrate

Art. ID TRC-A281810-10G
Unit 10 g
Deliverydetails Dry ice shipment /not restricted

Description

Category: Carbohydrates and Derivatives, Heterocycles, Metabolites and Impurities, Nucleotides, Bases and Related Reagent, Phosphorylating and Phosphitylating Agents /// Appearance: White to Off-White Solid /// Application Notes: A metabolite of Adenosine (A280400), a multifunctional nucleoside triphosphate used in cells as a coenzyme of intracellular energy transfer. It transports chemical energy within cells for metabolism. /// References: Knowles, J., et al.: Annu. Rev. Biochem., 49, 877 (1980)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine 5'-Triphosphate Disodium Salt Hydrate	[34369-07-8]				