

Aliquat 336

Art. ID TRC-A535950-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2810 III E1

Description

Application Notes: Aliquat 336 is a reagent in the synthesis of Iso Promethazine (I822500, HCl), an impurity of Promethazine (P757000, HCl) which is an antihistaminic, antiemetic, and CNS depressant. /// References: Gasco, M.R., et al.: J. Pharma. Sci., 68, 612 (1979)| Edge, W., et al.: J. Pharm. Pharmacol., 5, 279 (1953)| Huang, et al.: J. Pharm. Sci., 59, 772 (1970)| Shearer, C.M., et al.: Anal. Profiles Drug Subs., 5, 429 (1976)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tricaprylylmethylammonium chloride	[63393-96-4]				