

Allopurinol

Art. ID TRC-A547300-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2811 III E1

Description

Category: Heterocycles, Nucleotides, Bases and Related Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White to Off-White Solid /// Application Notes: Allopurinol is a xanthine oxidase inhibitor, decreases uric acid production. Used in treatment of hyperuricemia and chronic gout. Antiuric lithic. /// References: Benezra, S.A., et al.: Anal. Profiles Drug Subs., 7, 1, (1978)| Dalbeth, N., et al.: Rheumatology, 44, 1090 (2005)| Dinarello, C., et al.: J. Exp. Med., 201, 1355 (2005)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Allopurinol	[315-30-0]				