

## **D-allo-Threonine**

Art. ID TRC-A548540-5G  
Unit 5 g  
Deliverydetails No Dangerous Good

### **Description**

Category: Amines, Amino Acids and Derivatives, Chiral Reagents /// Application Notes: D-allo-Threonine is used to prepare serine and threonine beta-lactones as inhibitors of hepatitis A virus 3C cysteine proteinase. It is also used to synthesize the aglycon of ramoplanin A2 and ramoplanose with antibacterial activities. /// References: Lall, M., et al.: J. Org. Chem., 67, 1536 (2002) | Jiang, W., et al.: J. Am. Chem. Soc., 125, 1877 (2003)

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | D-allo-Threonine  | [24830-94-2] |                     |      |        |        |