

Aluminum Sulfate Hydrate

Art. ID TRC-A575735-500G
Unit 500 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Miscellaneous Reagents /// Application Notes: Aluminum Sulfate Hydrate is used in the synthesis of aluminum diethylphosphinate. Also used in the preparation of nanowires through a molten salt method. /// References: Yang, L. et al.: Chin. Chem. Lett., 22, 385 (2011)| Huo, K. et al.: J. Nanosci. Nanotechnol., 10, 4792 (2010)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminum sulfate hydrat e	[17927-65-0]				