

Aminoacetonitrile Hydrochloride

Art. ID TRC-A580055-100G
Unit 100 g
Deliverydetails No Dangerous Good

Description

Category: Miscellaneous Reagents /// Appearance: White to Off-White Solid /// Application Notes: Aminoacetonitrile is a useful synthetic intermediate. It is a reagent used to synthesize dipeptide nitriles as reversible and potent cathepsin S inhibitors. It can also be used to prepare substituted cyclic ureas as possible HIV protease inhibitors. /// References: Ward, Y., et al.: J. Med. Chem., 45, 5471 (2002)| Wilkerson, W., et al.: J. Med. Chem., 40, 4079 (1997)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aminoacetonitrile Hydrochloride	[6011-14-9]				