

Tin-base Lead-free Solder, Cast, Certified Reference Material (CRM), disc 38 mm x 13 mm

Art. ID MBH-74XGE2A
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,079 | % | | |
| | Aluminium (Al) | [7429-90-5] | 0,068 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,0086 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,713 | % | | |
| | Germanium (Ge) | [7440-56-4] | 0,479 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,031 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,0467 | % | | |