

Aminoguanidine Nitrate

Art. ID TRC-A609720-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1479 II E2

Description

Category: Miscellaneous Reagents /// Appearance: Off-White Solid /// Application Notes: Aminoguanidine Nitrate is used in the synthesis of inhibitors of cathepsin L used in tumor suppression and treatment. Also used in the synthesis of melanocortin-4 receptor agonists as potential antiobesity agents. /// References: Tian, X. et al.: J. Med. Chem., 51, 6055 (2008)| Song, J. et al.: Bioorg. Med. Chem. Lett., 23, 2801 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-Aminoguanidinium nitrate	[10308-82-4]				