

Antioxidant 24

Art. ID TRC-A697505-250G
Unit 250 g
Deliverydetails No Dangerous Good

Description

Appearance: White Solid /// Application Notes: Antioxidant 24 is used in the synthesis of composites to support polymer films, such as polypropylene and polyethylene film materials. /// References: del Teso Sanchez, K. et al.: Polym. Degrad. Stab., 113, 32 (2015)| Pantani, R. et al.: Eur. Polym. J., 49, 3471 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ultranox 626	[26741-53-7]				