

D-Arabinitol

Art. ID TRC-A735600-50G
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Carbohydrates and Derivatives, Diagnostics /// Application Notes: D-Arabinitol, can be used as a diagnostic tool for the detection of fungal infections in neutropenic patients. /// References: Salonen, J.H., et al.: Eur. J. Clin. Microbiol. Infect. Dis., 20, 179 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D(+)-Arabit	[488-82-4]				