

## **Arachidic Acid**

Art. ID TRC-A764990-1G  
Unit 1 g  
Deliverydetails No Dangerous Good /not restricted

### Description

Category: Aliphatics, Glycerols, Fatty Acid Derivatives and Lipids /// Appearance: White to Off-White Solid ///  
Application Notes: Arachidic acid is a saturated fatty acid and a minor constituent of peanut oil, coconut oil and corn oil. Arachidic acid is used for the production of detergents, photographic materials and lubricants. /// References: Kumar, S. R., et al.: J Chemtracks., 13, 277 (2011)| Poyarkova, T.N., et al: J. App. Chem. 85, 1259 (2012)|

| Text/Information | Analyte/Parameter                  | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------|------------|---------------------|------|--------|--------|
|                  | Arachidic acid, No. of<br>C: C20:0 | [506-30-9] |                     |      |        |        |