

## **Z-L-Arginine**

Art. ID TRC-A769400-25G  
Unit 25 g  
Deliverydetails No Dangerous Good /not restricted

### Description

Category: Amino Acids and Derivatives, Miscellaneous Reagents /// Appearance: White Solid /// Application Notes: Z-L-Arginine is used in the synthesis of renin inhibitors containing phosphorous acid derivatives at the scissile bond. Also used in the synthesis of a peptide deformylase inhibitor BB-3497. /// References: Allen, M. et al.: J. Med. Chem., 32, 1652 (1989)| Davies, S. et al.: Bioorg. Med. Chem. Lett., 13, 2715 (2003)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Z-Arg-OH	[1234-35-1]				