

Ascorbic Acid 2-sulfate Dipotassium Salt Dihydrate

Art. ID TRC-A787018-250MG
Unit 250 mg
Deliverydetails No Dangerous Good

Description

Appearance: White Solid /// Application Notes: Ascorbic Acid 2-sulfate Dipotassium Salt Dihydrate is an intermediate in the synthesis of L-threo-Hex-2-enaric acid 1,?4-Lactone which is an impurity of L-Ascorbic Acid (A786990) which is a physiological antioxidant. /// References: Al-Meshal, I.A., et al.: Anal. Profiles Drug Subs., 11, 45 (1982), Levine, M., et al.: N. Engl. J. Med., 314, 892 (1986), Prust, C., et al.: Nature Biotech., 23, 195 (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ascorbic Acid 2-sulfate Dipotassium Salt Dihyd rate					