

D,L-Azetidine-2-carboxylic Acid-d4

Art. ID TRC-A812492-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Application Notes: D,L-Azetidine-2-carboxylic Acid-d4 is the isotope labelled analog of D,L-Azetidine-2-carboxylic Acid (A812490), a four-membered ring analog of L-Proline and a useful intermediate in the synthesis of polypeptides. /// References: Barber, M., et al.: Int. J. Pept. Protein Res., 14(3), 247 (1979)| Tsai., F.H., et al.: Biopolymers, 30(11-12), 1039 (1990)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D,L-Azetidine-2-carboxy lic Acid-d4					