

Azobisisobutyronitrile

Art. ID TRC-A932550-100G

Unit 100 g

Description

Category: Miscellaneous Reagents /// Appearance: White Solid /// Application Notes: Azobisisobutyronitrile is used in the optimization of benzodiazepinones as selective inhibitors of apoptosis protein (XIAP). Also used in the preparation of novel stimuli-responsive hydrogels based on starch. /// References: Kester, R. et al.: J. Med. Chem., 56, 7788 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Azobisisobutyronitrile					