

## **N-Acetyl-D-alanine-2,3,3,3-d4**

Art. ID TRC-A167966-100MG  
Unit 100 mg  
Deliverydetails No Dangerous Good

### Description

Application Notes: N-Acetyl-D-alanine-2,3,3,3-d4 is the deuterium labelled form of N-Acetyl-DL-alanine (A167965), which is a derivative of DL-Alanine (A481400) and is used as a reagent in the synthesis of potent new agents for the treatment of epilepsy. /// References: Conley, J.D., et al.: J. Med. Chem., 30, 567 (1987)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-D-alanine-2,3, 3,3-d4					