

## **D-Alanine Isopropyl Ester Hydrochloride**

Art. ID TRC-A510820-25G  
Unit 25 g  
Deliverydetails No Dangerous Good

### **Description**

Appearance: White to Off-White Solid /// Application Notes: D-Alanine Isopropyl Ester Hydrochloride is a chiral reagent used in the synthesis of nucleotide prodrugs such as (PSI-7977) for the treatment of hepatitis C. Also used in the preparation of synthetic sweeteners. /// References: Sofia, M. et al.: J. Med. Chem., 53, 7202 (2010)| Mazur, R. et al.: J. Med. Chem., 16, 1284 (1973)|

| Text/Information | Analyte/Parameter                       | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------------|--------------|---------------------|------|--------|--------|
|                  | D-Alanine Isopropyl Ester Hydrochloride | [39613-92-8] |                     |      |        |        |